

# FITOVERS BARRIER TECHNOLOGIES

## Fitovers

- Perfect fit over existing personal Rx eyewear
- Goggle design protects from fluid splashes

<b>Protection:</b>	0.75mm Pb Equivalent Corning Med-X Lenses
<b>Side Shields:</b>	Yes – 0.50mm Pb Equivalent
<b>Lens Type:</b>	Plano only
<b>Frame Colour:</b>	Black, Tortoise
<b>Weight:</b>	88g

[Request a quote](#)

