

GLOBAL FIT GLASSES

RUDY PROJECT

Rudy Project Contender

- Lightweight titanium and carbon frame
- Ergonose system for a superior fit

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Black

Weight: 85g



[Request a quote](#)

Rudy Project Sensor

- Ultralight hybrid titanium and kynetium frame
- Ergonose system for an extraordinary fit

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano, single vision, bifocal, progressive

Frame Colour: Carbon silver

Weight: 54g



[Request a quote](#)

Rudy Project Voltage

- Ultralight titanium frame
- Ergonose system for maximum comfort

Protection: 0.75mm Pb Equivalent
Corning Med-X Glass

Side Shields: None

Lens Type: Plano only

Frame Colour: Carbon grey

Weight: 59g



[Request a quote](#)

GLOBAL FIT GLASSES

OAKLEY

Oakley Crosslink

- Sleek and modern look with a comfortable fit
- Constructed of lightweight, durable O-Matter™ material

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano, Single Vision,
BiFocal, Progressive

Frame Colour: Grey Smoke, Satin Black

Weight: 58g



[Request a quote](#)

Oakley Straightlink

- O-Matter™ stress-resistant frame for durability and all-day comfort
- Unobtainium® earsocks and nosebombs increase grip with perspiration for a no slip grip without compromising on style

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano only

Frame Colour: Matte Black

Weight: 66g



[Request a quote](#)

Oakley TwoFace

- Unique pairing of lightweight O-Matter™ upper and high strength aluminium lower rims
- Lens positioning offers superior comfort, clarity and razor sharp vision at every angle

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: None

Lens Type: Plano, Single Vision,
BiFocal, Progressive

Frame Colour: Polished Black with
Steel Grey, Triple Black

Weight: 64g



[Request a quote](#)

GLOBAL FIT GLASSES BARRIER TECHNOLOGIES

Blade II

- Adjustable nose piece for a comfortable fit
- Integrated side shields for enhanced lateral protection

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano, Single Vision

Frame Colour: Black

Weight: 64g



[Request a quote](#)

MX-30

- Super soft rubber nose piece
- Global Fit

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano, Single Vision, BiFocal, Progressive

Frame Colour: Black, Blue, Crimson, Silver

Weight: 69g



[Request a quote](#)