

NON-PRESCRIPTION LEAD GLASSES

RUDY PROJECT

Rudy Project Contender

- Lightweight titanium and carbon frame
- Ergonose system for a superior fit

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Black

Weight: 85g



[Request a quote](#)

Rudy Project Sensor

- Ultralight hybrid titanium and kynetium frame
- Ergonose system for an extraordinary fit

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano, single vision, bifocal, progressive

Frame Colour: Carbon silver

Weight: 54g



[Request a quote](#)

Rudy Project Voltage

- Ultralight titanium frame
- Ergonose system for maximum comfort

Protection: 0.75mm Pb Equivalent
Corning Med-X Glass

Side Shields: None

Lens Type: Plano only

Frame Colour: Carbon grey

Weight: 59g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

OAKLEY

Oakley Crosslink

- Sleek and modern look with a comfortable fit
- Constructed of lightweight, durable O-Matter™ material

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano, Single Vision, BiFocal, Progressive

Frame Colour: Grey Smoke, Satin Black

Weight: 58g



[Request a quote](#)

Oakley Fives Squared

- O-Matter™ stress-resistant frame
- Comfortable and durable for all day comfort and protection

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: None

Lens Type: Plano only

Frame Colour: Polished Black

Weight: 58g



[Request a quote](#)

Oakley Fuel Cell

- Unobtanium™ sure-grip inset zones provide durability and performance
- Three-point fit holds lenses in precise optical alignment

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Polished Black, Multi-Camo

Weight: 82g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

OAKLEY

Oakley Gascan

- O-Matter™ stress-resistant frame
- Lightweight and durable for all-day comfort and protection

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Matte Black, Polished Black

Weight: 74g



[Request a quote](#)

Oakley Holbrook

- O-Matter™ stress-resistant frame
- Lightweight and durable for all-day comfort and protection

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: None

Lens Type: Plano, Single Vision, BiFocal, Progressive

Frame Colour: Satin Black, Polished Clear

Weight: 65g



[Request a quote](#)

Oakley Jupiter Squared

- Unique Oakley riveting on temple arm
- O-Matter™ frame ensures comfort during long procedures

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano, Single Vision, BiFocal, Progressive

Frame Colour: Woodgrain, Polished Black

Weight: 65g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

OAKLEY

Oakley Straightlink

- O-Matter™ stress-resistant frame for durability and all-day comfort
- Unobtainium® earsocks and nosebombs increase grip with perspiration for a no slip grip without compromising on style

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano only

Frame Colour: Matte Black

Weight: 66g



[Request a quote](#)

Oakley TwoFace

- Unique pairing of lightweight O-Matter™ upper and high strength aluminium lower rims
- Lens positioning offers superior comfort, clarity and razor sharp vision at every angle

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: None

Lens Type: Plano, Single Vision, BiFocal, Progressive

Frame Colour: Polished Black with Steel Grey, Triple Black

Weight: 64g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

COSTA DEL MAR

Costa del Mar Blackfin

- Durable lightweight TR-90 nylon frame
- Hydrolite hypo-allergenic rubber prevents slippage

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Matte Black

Weight: 73g



[Request a quote](#)

Costa del Mar Brine

- Durable lightweight TR-90 nylon frame
- Vented ports provide maximum airflow and prevent fogging

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: None

Lens Type: Plano, Single Vision

Frame Colour: Black, Tortoise

Weight: 56g



[Request a quote](#)

Costa del Mar Cape

- Hydrolite hypo-allergenic rubber prevents slippage
- Vented ports promote maximum airflow and prevents fogging

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Matte Black

Weight: 92g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

COSTA DEL MAR

Costa del Mar Fantail

- Modified DBL (distance between lenses) eliminates frame resting on face
- Hydrolite hypo-allergenic rubber prevents slippage

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Black

Weight: 69g



[Request a quote](#)

Costa del Mar Permit

- Modified DBL (distance between lenses) eliminates frame resting on face
- Vented ports promote maximum airflow and prevents fogging

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Matte Black

Weight: 69g



[Request a quote](#)

Costa del Mar Tuna Alley

- Modified DBL (distance between lenses) eliminates frame resting on face
- Vented ports promote maximum airflow and prevents fogging

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Matte Black, Tortoise

Weight: 69g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

BARRIER TECHNOLOGIES

Astro II

- Integrated peripheral splash shields for superior splatter protection
- Rubberised nose pads for comfort

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano, Single Vision, BiFocal, Progressive

Frame Colour: Blue, Pink

Weight: 67g



[Request a quote](#)

Atomic

- Vented frame style prevents fogging
- Rubberised anti-slip frame

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: Yes – 0.50m Pb Equivalent

Lens Type: Plano, Single Vision, BiFocal, Progressive

Frame Colour: Black, Gunmetal

Weight: 60g



[Request a quote](#)

Barrier Lite

- Lightweight polycarbonate frame
- Wrap around style

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano, Single Vision

Frame Colour: Black, Silver

Weight: 60g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

BARRIER TECHNOLOGIES

Barrier Brazen

- Universally fitting frame for all face sizes
- Ventilated nose bridge for additional comfort

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano, Single Vision,
BiFocal, Progressive

Frame Colour: Black, Tortoise

Weight: 86g



[Request a quote](#)

Blade II

- Adjustable nose piece for a comfortable fit
- Integrated side shields for enhanced lateral protection

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: Yes – 0.50mm Pb
Equivalent

Lens Type: Plano, Single Vision

Frame Colour: Black

Weight: 64g



[Request a quote](#)

Circuit XL

- Anti-slip rubber nose pads and adjustable temples
- Fits most face types

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb
Equivalent

Lens Type: Plano, Single Vision,
BiFocal, Progressive

Frame Colour: Black, Metallic Blue, Red

Weight: 65g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

BARRIER TECHNOLOGIES

Fitovers

- Perfect fit over existing personal Rx eyewear
- Goggle design protects from fluid splashes

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Black, Tortoise

Weight: 88g



[Request a quote](#)

Grid II

- Wrap frame provides excellent lateral protection
- Ultimate in comfort and attenuation

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano, Single Vision, BiFocal, Progressive

Frame Colour: Gloss Black, Gunmetal Grey, Navy

Weight: 70g



[Request a quote](#)

MX-30

- Super soft rubber nose piece
- Global Fit

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano, Single Vision, BiFocal, Progressive

Frame Colour: Black, Blue, Crimson, Silver

Weight: 69g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES BARRIER TECHNOLOGIES

Slim-T

- Precision machined ultralight frame
- Dual hinges for ultimate comfort and flex

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: Yes – 0.50mm Pb Equivalent

Lens Type: Plano only

Frame Colour: Gunmetal Grey

Weight: 79g



[Request a quote](#)

T-Zone

- Fashionable wrap around style
- Available in 3 sizes – Small, Medium, Large


Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: Yes – 0.50m Pb Equivalent

Lens Type: Plano, Single Vision,
BiFocal, Progressive

Frame Colour: Black, Blue, Red, Silver

Weight: 60g (SM) / 66g (MD) /
80g (LG)



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

RAY-BAN

Ray-Ban Balorama

- Wrap-around style
- Great coverage

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano only

Frame Colour: Black

Weight: 74g



[Request a quote](#)

Ray-Ban Jackie Ohh

- Retro classic shape – Hollywood glamour!
- Great coverage


Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: None

Lens Type: Plano only

Frame Colour: Black

Weight: 69g



[Request a quote](#)

Ray-Ban Justin

- Inspired by the iconic Wayfarer
- Slightly larger rectangular lenses

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: None

Lens Type: Plano only

Frame Colour: Black

Weight: 65g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

RAY-BAN

Ray-Ban New Wayfarer

- Update on the classic Wayfarer
- Lightweight protection

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano, Single Vision, BiFocal,
Progressive

Frame Colour: Black

Weight: 59g



136mm



18mm



55mm



41mm



145mm

[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

WILEY X

Wiley X P-17

- Double injected rubber nose piece for comfort
- Lightweight design

Protection: 0.75mm Pb Equivalent
Corning Med-X Lenses

Side Shields: None

Lens Type: Plano only

Frame Colour: Gloss Black

Weight: 56g



[Request a quote](#)

Wiley X Gravity

- Flexible half rubber temples
- Wrap-around style


Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: Yes – 0.50mm Pb
Equivalent

Lens Type: Plano only

Frame Colour: Matte Black

Weight: 76g



[Request a quote](#)

Wiley X Brick

- Wrap-around style
- Lightweight frame

Protection: 0.75mm Pb Equivalent
Corning Med-X lenses

Side Shields: None

Lens Type: Plano only

Frame Colour: Matte Black, Crystal
Metallic

Weight: 56g



[Request a quote](#)

NON-PRESCRIPTION LEAD GLASSES

WILEY X

Wiley X Enzo

- Wrap-around frame
- Larger lens for greater coverage

Protection:	0.75mm Pb Equivalent Corning Med-X Lenses
Side Shields:	Yes – 0.50mm Pb Equivalent
Lens Type:	Plano only
Frame Colour:	Matte Black
Weight:	66g



135mm



15mm



64mm



41mm



125mm

[Request a quote](#)