

# PRESCRIPTION LEAD GLASSES

## RUDY PROJECT

### Rudy Project Sensor

- Ultralight hybrid titanium and kynetium frame
- Ergonose system for an extraordinary fit

<b>Protection:</b>	0.75mm Pb Equivalent Corning Med-X Lenses
<b>Side Shields:</b>	None
<b>Lens Type:</b>	Plano, single vision, bifocal, progressive
<b>Frame Colour:</b>	Carbon silver
<b>Weight:</b>	54g



138mm



14mm



64mm



37mm



135mm

[Request a quote](#)

# PRESCRIPTION LEAD GLASSES

## OAKLEY

### Oakley Crosslink

- Sleek and modern look with a comfortable fit
- Constructed of lightweight, durable O-Matter™ material

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** None

**Lens Type:** Plano, Single Vision,  
BiFocal, Progressive

**Frame Colour:** Grey Smoke, Satin Black

**Weight:** 58g



[Request a quote](#)

### Oakley Holbrook

- O-Matter™ stress-resistant frame
- Lightweight and durable for all-day comfort and protection

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** None

**Lens Type:** Plano, Single Vision,  
BiFocal, Progressive

**Frame Colour:** Satin Black, Polished Clear

**Weight:** 65g



[Request a quote](#)

### Oakley Jupiter Squared

- Unique Oakley riveting on temple arm
- O-Matter™ frame ensures comfort during long procedures

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** None

**Lens Type:** Plano, Single Vision,  
BiFocal, Progressive

**Frame Colour:** Woodgrain, Polished Black

**Weight:** 65g



[Request a quote](#)

# PRESCRIPTION LEAD GLASSES OAKLEY

## Oakley TwoFace

- Unique pairing of lightweight O-Matter™ upper and high strength aluminium lower rims
- Lens positioning offers superior comfort, clarity and razor sharp vision at every angle



<b>Protection:</b>	0.75mm Pb Equivalent Corning Med-X lenses
<b>Side Shields:</b>	None
<b>Lens Type:</b>	Plano, Single Vision, BiFocal, Progressive
<b>Frame Colour:</b>	Polished Black with Steel Grey, Triple Black
<b>Weight:</b>	64g



140mm



16mm



60mm



41mm



135mm

[Request a quote](#)

# PRESCRIPTION LEAD GLASSES

## COSTA DEL MAR

### Costa del Mar Brine

- Durable lightweight TR-90 nylon frame
- Vented ports provide maximum airflow and prevent fogging

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X lenses

**Side Shields:** None

**Lens Type:** Plano, Single Vision

**Frame Colour:** Black, Tortoise

**Weight:** 56g



127mm



18mm



59mm



36mm



118mm

[Request a quote](#)

# PRESCRIPTION LEAD GLASSES

## RAY-BAN

### Ray-Ban New Wayfarer

- Update on the classic Wayfarer
- Lightweight protection

<b>Protection:</b>	0.75mm Pb Equivalent Corning Med-X Lenses
<b>Side Shields:</b>	None
<b>Lens Type:</b>	Plano, Single Vision, BiFocal, Progressive
<b>Frame Colour:</b>	Black
<b>Weight:</b>	59g



136mm



18mm



55mm



41mm



145mm

[Request a quote](#)

# PRESCRIPTION LEAD GLASSES BARRIER TECHNOLOGIES

## Astro II

- Integrated peripheral splash shields for superior splatter protection
- Rubberised nose pads for comfort

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** Yes – 0.50mm Pb Equivalent

**Lens Type:** Plano, Single Vision, BiFocal, Progressive

**Frame Colour:** Blue, Pink

**Weight:** 67g



[Request a quote](#)

## Atomic

- Vented frame style prevents fogging
- Rubberised anti-slip frame

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X lenses

**Side Shields:** Yes – 0.50m Pb Equivalent

**Lens Type:** Plano, Single Vision, BiFocal, Progressive

**Frame Colour:** Black, Gunmetal

**Weight:** 60g



[Request a quote](#)

## Barrier Lite

- Lightweight polycarbonate frame
- Wrap around style

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** Yes – 0.50mm Pb Equivalent

**Lens Type:** Plano, Single Vision

**Frame Colour:** Black, Silver

**Weight:** 60g



[Request a quote](#)

# PRESCRIPTION LEAD GLASSES BARRIER TECHNOLOGIES

## Barrier Brazen

- Universally fitting frame for all face sizes
- Ventilated nose bridge for additional comfort

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** Yes – 0.50mm Pb Equivalent

**Lens Type:** Plano, Single Vision,  
BiFocal, Progressive

**Frame Colour:** Black, Tortoise

**Weight:** 86g



[Request a quote](#)

## Blade II

- Adjustable nose piece for a comfortable fit
- Integrated side shields for enhanced lateral protection

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X lenses

**Side Shields:** Yes – 0.50mm Pb  
Equivalent

**Lens Type:** Plano, Single Vision

**Frame Colour:** Black

**Weight:** 64g



[Request a quote](#)

## Circuit XL

- Anti-slip rubber nose pads and adjustable temples
- Fits most face types

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** Yes – 0.50mm Pb  
Equivalent

**Lens Type:** Plano, Single Vision,  
BiFocal, Progressive

**Frame Colour:** Black, Metallic Blue, Red

**Weight:** 65g



[Request a quote](#)

# PRESCRIPTION LEAD GLASSES BARRIER TECHNOLOGIES

## Grid II

- Wrap frame provides excellent lateral protection
- Ultimate in comfort and attenuation

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** Yes – 0.50mm Pb Equivalent

**Lens Type:** Plano, Single Vision, BiFocal, Progressive

**Frame Colour:** Gloss Black, Gunmetal Grey, Navy

**Weight:** 70g



[Request a quote](#)

## MX-30

- Super soft rubber nose piece
- Global Fit

**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** Yes – 0.50mm Pb Equivalent

**Lens Type:** Plano, Single Vision, BiFocal, Progressive

**Frame Colour:** Black, Blue, Crimson, Silver

**Weight:** 69g



[Request a quote](#)

## T-Zone

- Fashionable wrap around style
- Available in 3 sizes – Small, Medium, Large


**Protection:** 0.75mm Pb Equivalent  
Corning Med-X Lenses

**Side Shields:** Yes – 0.50mm Pb Equivalent

**Lens Type:** Plano, Single Vision, BiFocal, Progressive

**Frame Colour:** Black, Blue, Red, Silver

**Weight:** 60g (SM) / 66g (MD) / 80g (LG)



[Request a quote](#)